ALE-HOP

User manual

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Manual d'usuari

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Mode d'emploi

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Manuale d'uso

Priručnik s uputama

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Gross weight: 123 grams

Power: 2x AA batteries (not included)

1.1 Concept of Low-energy Infrared Hair Growth

Low-energy infrared is a monochrome wavelength capable of reaching its climax to enable it to infiltrate into substances without any heat effect or causing any damage or side effects.

There are two kinds of infrared light: high-energy infrared, which can penetrate subcutaneous tissue and generate heat, and low-energy infrared, which has no heat nor function of light or power.

Modern science research shows that infrared can penetrate superficial tissue to raise the activity of ATP, the substance responsible for transmitting energy among cells. Therefore, low-energy infrared helps activate cells, influence tissue replacement, accelerate reorganization, and assist in stopping inflammation. The infrared hair growth comb uses light to clean your hair roots and makes hair shine.

1.2 Concept of O₃ dandruff removing

Dandruff comes from fungus generated massively on the surface of hair when the hair does not dry enough after being washed. These fungi absorb sebum and quickly expand. The sebum starts to produce many secretions and increases the alternate rate of epithelial tissue. It is the repeated vicious

cycle that causes dandruff.

O₃ can kill fungus in humid environments. Use the massage comb while your hair is a bit wet. After three weeks of using the massage comb, you will notice less dandruff in your hair.

1.3 Concept of shock massage

Besides advanced hair growth and care by infrared, the hair comb has a shock massage function. It refers to massage on the scalp by an external force, which helps blood circulate, relax, and relieve fatigue. Frequent use can prevent headaches or intermittent headaches.

Scientific theory proves that the liquid of the soft tissue in the human body has a similar frequency to radiated magnetic waves and infrared. One cell subsists with health transmittance by relying on cells surrounding it. Each cell completes its process with an individual wave motion system influenced by a certain quantity of low-energy infrared lights that transfer messages through interconnected particle areas.

2. Infrared hair comb function

Infrared Hair Growth Comb is an instrument made for hair and scalp care. The laser has a wavelength of 650 nm that can activate atrophic hair follicles, balance biological energy, and promote metabolism. The scalp and follicle tissue can fully absorb nutrients in a normal metabolism. Moreover, infrared helps blood circulation and fastens the process of nutrients reaching the follicle tissue. Hair growth can be stimulated by accelerating the cycle in its early, middle, and end stages.

HOW TO USE

1. Dry your hair to an 80% of extent.
2. Open the battery cover to plug in 2x AA batteries and close off the battery compartment.
3. Press the ON/OFF button. Two infrared laser lights will come out from the comb.
4. Place the comb on the head and start combing from the forehead hairline to the back of the head for about 10 minutes.

MESSAGE COMB FUNCTIONS

1. Stimulates the hair follicle

Low-energy infrared can effectively stimulate hair follicles and activate 83% of telogen hair follicles into anagen follicles by continuous radiation treatment, thus preventing hair loss and promoting hair growth.

2. Accelerates blood circulation and promotes metabolism

Low-energy infrared can increase blood circulation by 20-30% in the front head and 40-60% in the head top. The laser has been tested to help boost blood circulation and send more oxygen and nutrients to the hair follicles to prevent hair loss.

3. Controls oil

When the hair secretes excessive grease is mainly caused by sebaceous glands that cause the hair to become oiled even after one day after being cleaned and could end up in hair loss. Soft laser can make sebaceous glands smaller and decrease or control oil to improve hair loss.

4. Promotes hair health

Frequent use of low-energy infrared light in radiating hair follicles can strengthen the hair's tightness, proximity, and flexibility

and accelerate hair growth.

WARNINGS

When using electric products, basic precautions should always be followed, including the following:

1. Read the instructions before the first use of the product.
2. O₃ produces a smelly gas. Please be assured that the slight amount of O₃ from the product will not cause any harm to you.
3. Never stare directly or point the O₃ export light at someone else's face for a long time. If the light turns dark, please replace or charge the battery.
4. If you do not plan on using the massage comb for a long time, please extract the battery from the comb so the battery power is not consumed.
5. Do not use the device around inflammable liquids such as benzene or fire extinguishers. There could be a risk of explosion or fire.
6. Do not use the massage comb with wet hands. It could cause a power short circuit.
7. This product is not suitable for children or babies. Please consider this when using the product.
8. Do not use diluents, benzene, or other solvents to clean the product. They could damage the product.

BOX CONTENT

- 1x Infrared electric massage comb
- 1x Instruction manual

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Manufactured by | **CLAVE DÉNIA, S.A.U.**