
## ALE-HOP

**User manual** | Manual de usuario | Manual  
d'usuari | Manual do usuário | Mode  
d'emploi | Benutzerhandbuch | Manuale  
d'uso | Priručnik s uputama | للتحويل بين الأوضاع

# RECHARGEABLE ENERGY MASSAGE PEN

## MASSAGE HEADS

Node point massage

Traditional  
Acupuncture

Facial Beauty  
Acupuncture

**Node point massage:** Node therapy for partial body pains

**Traditional acupuncture:** Partial therapy for quickly alleviating pain.

**Facial beauty:** Superficial therapy for health care and facial beauty.

## CHARGING THE DEVICE

Plug the charging cable into the device and the other end into a power source. The red charging indicator light will turn on. Please charge the device for 2 hours before its first use.

Once the device is fully charged, the red charging indicator light will turn off.

## HOW TO USE

Turn on the device by pressing the power button. Press the + button to start the massage pen. The number 1 will appear on the display. To increase the intensity, press the + button (max level 9; To decrease the intensity, press the - button (min level 0). Long press the power button for 3 seconds to turn off the device.

### Laser open method

Turn on the device. The numbers 1 – 9 will appear on the display. Press the power button once or twice to open the laser and press three times to close the laser.

Start with the device on a medium-low setting (1 or 2). If you do not feel anything, slide the massage head around the general area until you feel the pulsing sensation. If you do not feel the pulsing sensation, press the + button to a higher setting.

Continue sliding the massage head around until you feel a strong pulse. When the pulse is strong, you will have found an acupoint that is being stimulated. You may hold the massage pen directly on the point or gently massage the acupoint and the associated meridian. It is possible to feel the pulsing on another point along the meridian or coursing through the entire length of the meridian. Locations of points and meridians are at the end of the instructions.

### ACUPOINTER AS A TENS UNIT

Using the massage pen to trigger muscle contraction is a great way to exercise a hard-to-isolate muscle, to help loosen a frozen muscle, or to prevent or reverse the atrophy of immobilized muscles when you have a broken bone in a cast. (Consult your doctor before using the device directly on a wound or broken bone.)

It takes a little practice to find the muscle's trigger points, so the device must be pushed a little deeper into the surface of the skin. It is also important to relax the muscle you are trying to trigger as tensing may interfere with the device.

When you find the right spot, the muscle will "jump" and "bounce" vigorously until the tip moves from that spot. A slight movement will stop the effect, so try to find a comfortable position to stand for several minutes (fifteen to twenty minutes is optimal). Now relax the limb and let the massage pen perform the isolated muscle workout.

### MESSAGE PEN AS A PAIN RELIEVER

Apply aloe vera gel or any moisturizer to the affected area and press the massage head against the skin, making firm contact. Hold the device in one spot or move it until you feel a strong pulsation. Continue for several minutes.

Body hair can interfere with the device making good contact with the skin. Using more aloe vera gel or any moisturizer will prevent the electrical arcing that causes a stinging sensation. Use a water-based analgesic balm to use the massage pen if you have muscle pain.

## **JOLT THE JOINTS**

When experimenting with the massage pen, a great place to start is at the joints: elbows, knees, ankles, and wrists. The joints are one of the parts where external pathogens linger in the body and stay trapped, causing blockages of Qi. Use low-voltage therapy in these areas to help the body in purging pathogens.

## **MESSAGE PEN AS A BEAUTY TREATMENT**

Set the massage pen on a lower setting (2 to 4) and apply aloe vera gel or any cleansing or moisturizing product to the face, as long as the face remains wet during the treatment. Gently, but still making contact, move the massage pen around the face but avoid contact with eyes or mucous membranes. When you feel pulsing, leave it on the spot for a few minutes and let the device exercise, tone, and tighten those facial muscles that tend to sag with age.

Concentrate on problematic spots such as bags or puffiness around the eyes and cheeks. You'll be amazed at how well it works after a few treatments!


## **BENEFITS OF THE MESSAGE PEN**

1. The massage pen can be used as a portable hand-held "TENS" unit to automatically trigger rapid muscle contraction when muscles need to be isolated or when the muscle is immobilized in a cast or splint. The pulsed contraction effectively works the muscle for you to prevent atrophy.
2. The pulsations from the massage pen can relieve pain by blocking pain signals as they travel from the nervous system to the brain. They also have a relaxing effect on tense muscles as a massage.
3. As a beauty aid for stimulating the muscles of the face and body as well as tightening and rejuvenating the skin and subcutaneous layers to reduce wrinkles and sagging flesh associated with aging.
4. The massage pen helps locate acupoints on the body and stimulates them to unblock Qi flow along a meridian or to focus Qi on an acupoint. The mild electrical pulses act in the same way as acupuncture needles or acupoint massage.

## **ACUPOINTS AND MERIDIANS**


The following images show the locations of all the major meridians and acupoints along with a list of common ailments for which they can be

used. Note the golden coloured dots are the Golden Acupoints, the most potent and effective when stimulated with the massage


### HEART MERIDIAN

Excessive sweating, heart palpitations, fatigue, insomnia, poor circulation in the extremities, dizziness, mouth sores and nightmares


### BLADDER MERIDIAN

Joint pain, lower back pain, constipation, fatigue, rapid pulse, tinnitus, decreased sex drive, reproductive issues


### SMALL INTESTINE MERIDIAN

Fever, chronic cough, shoulder and neck pain, sore throat, TMJ syndrome, swollen glands, poor circulation in extremities


### GALL BLADDER MERIDIAN

Headaches, migraine, rib-side pain, swollen knees, leg weakness and numbness, chronic fatigue, anxiety, insomnia


### LIVER MERIDIAN

Mood disorders, headache, nausea, muscle atrophy, cramps, spasms


### LUNG MERIDIAN

Asthma, common cold, allergies, stress, acne, eczema, coughing, sore throat, dry skin


### KIDNEY MERIDIAN

Sexual dysfunction, impotence, tinnitus, oedema, premature greying, vertigo, respiratory ailments


### PERICARDIUM MERIDIAN

Anxiety attacks, nausea, morning sickness, motion sickness, depression, dry mouth


### TRIPLE BURNER MERIDIAN

Depression, mood swings, irritability, sore throat, earaches, tinnitus, headaches, neck and head pain.


### LARGE INTESTINE MERIDIAN

Fever, allergies, headaches, sneezes, congestion, colds, flu, itching, constipation, upper back pain

### STOMACH MERIDIAN

Acid reflux, indigestion, nausea, toothache, muscle spasms, diarrhoea, constipation

### SPLEEN MERIDIAN

Digestive disorders, dyspepsia, nausea, acid reflux, constipation diarrhoea

## **WARNINGS**

When using electric products, basic precautions should always be followed, including the following:

1. Read all the instructions before using the product.
2. Do not use the massage pen if you are pregnant, have a life-threatening disease or injury, or have a heart condition and a sustaining electronic medical device such as a pacemaker, pulse regulator, artificial heart or lungs.
3. Never use the massage pen in the shower, bath, pool, or high-humidity environments.
4. Do not use on or near mucous membranes or in any orifice of the body.
5. Keep away from children and pets. Do not use it on anyone who is unconscious or is unable to speak.
6. The message pen uses your hand as a ground pole to complete the circuit, so you may sometimes feel the voltage in the hand you are holding it with. That's because there are many acupoints in the hands.
7. Change your grip slightly or grip the device with a larger area of your hand to disperse the effect.
8. Do not disassemble, use, or try to repair the product if it is damaged.

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Manufactured by | **CLAVE DÉNIA, S.A.U.**